

OXIDIZING SOLIDS

TEST APPARATUS TEST FOR OXIDIZING SOLIDS

Model: BB-12

ver. 1.3, 2020




STANDARD REFERENCES:

- UN Recommendation on the Transport of Dangerous Goods.
Manual of Tests and Criteria,
Test O.1: Test for oxidizing solids.



SPECIFICATION:

Gloving wire length:	30 cm
Gloving wire diameter:	Ø0,6 mm
Wire heating power:	150W
Insulating plate :	150x150 mm
Glass plate:	quartz glass
Plate holder:	stainless steel
Glass funnel:	60°, 70mm
Power supply:	110 / 230 VAC
Remote switch:	on/off
Wire forming rig:	as per O.1 test



ACCESSORIES:

Spare glass plate:	2 pcs.
Heating wire:	100 m
Spoons:	1 pc.
Spatula:	1 pc.
Stop-watch:	1 pc.
Small brush:	1 pc.
Manuals:	User Manual Maintenance Guide

The information given in this document represents the state of engineering at the time of publishing. We reserve the right to make modifications to above specifications.